

LIFEPAK® CR2

cellular connected defibrillator

Sudden cardiac arrest (SCA) is a leading cause of death in Australia and New Zealand. Timely response by first responder teams and early intervention by bystanders remains a key factor driving favourable outcomes for victims of SCA. That's why public access defibrillators are so important!

Anyone can use a defibrillator and start CPR to try and save a life of someone in cardiac arrest. But AEDs are only effective if they are easily accessible and ready to be used in a sudden cardiac arrest emergency.



Benefits of a connected AED solution

Reduce risk



Are you prepared? Ensure your AED is ready to be used in an SCA emergency. Self-monitoring AED completes automatic self-tests and displays state of readiness in LIFELINKcentral™ AED program manager. Monitor AED's location if it has been moved, lost or stolen.

Used your device? Remote cellular connectivity enables real-time notifications if AED has been used or requires servicing.

Save time and money



Improve operational efficiency. Automate fleet management tasks by remotely managing expiry dates and replacement needs. Reduce risk of incorrect data entry when manually tracking AEDs.

Increase return on investment. Reduce human resources and manual labour time required to manage AED fleet.

Stay secure



All data is stored in Australia. Data is kept secure using a custom-built cloud infrastructure.

In case of an emergency:

- 1 Call 000
- 2 Begin CPR until medical treatment arrives
- 3 Use AED and follow verbal or visual prompts*

*An AED will only deliver a shock to a person with a shockable rhythm

Cellular and Wi-Fi connectivity available

Self-monitoring AED tracks readiness on one platform

Get real-time notifications when device is used or requires servicing

Automatically track expiry dates for AED accessories

Track geolocation of AEDs

Enable AED Locator feature to locate AED if it is moved outside of its set diameter

Remotely manage automatic software updates

Tailor LIFELINKcentral AED program manager to suit business needs

ALWAYS FOLLOW THE DIRECTIONS FOR USE.

The information presented is intended to demonstrate Stryker's product offerings. Users must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: LIFEPAK, LIFELINKcentral, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

CA4351_SP © 2024 Stryker



Manufactured by:
Physio-Control, Inc.
11811 Willows Road NE Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
stryker.com/emergencycare

Stryker Australia
8 Herbert Street, St Leonards NSW 2065
T: 61 2 9467 1000 F: 61 2 9467 1010

Stryker New Zealand
511 Mount Wellington Hwy,
Mt Wellington Auckland 1060
T: 64 9 573 1890 F: 64 9 573 1891